

2026年2月

お得意先様 各位

アズワン株式会社
営業本部

価格改定のご案内

拝啓 貴社益々ご清栄こととお慶び申し上げます。平素は格別なるご高配を賜り、厚く御礼申し上げます。

掲題の件、原材料コストや輸送費の上昇等により、誠に不本意ではございますが、
一部の商品につき、2026年3月1日付で価格改定を実施させて頂きたく、ご案内申し上げます。

弊社といたしましても、ユーザー様、販売店様へご迷惑の掛からぬよう企業努力で対応して参りましたが、
弊社の努力だけでは吸収しがたい状況となつてきております。

ユーザー様、販売店様にはご迷惑をお掛けいたしますが、諸般の事情をご賢察頂き、
ご理解賜りますようお願い申し上げます。

敬具

価格改定表

研究用総合機器カタログ2025
サンクアスト産業用研究機器カタログ2025
ナビス看護・医療用品総合カタログ2026
サニーフーズカタログ2024
研究設備カタログ2024
アズピュア クリーン環境・静電気対策カタログ2025
アズツール 生産現場向け副資材カタログ2021
サイエンスカタログ2024-2025
災害対策カタログ2025

2026年3月1日

 **アズワン株式会社**

品番	商品名	旧価格	新価格	研究用総合 機器2025 頁	サンクアスト 2025 頁	ナビス2026 頁	サニーフー 2024 頁	設備2024 頁	アズビエア 2025 頁	アズツール 2021 頁	サイエンスカ タログ 頁	災害対策カ タログ 頁
0-1506-01	殺菌ランプ(フィリップス社) TUV-15W	3,980	4,980			1033						
0-5167-01	鉗子立カバー 1号 BMC-8960	860	950			1186						
0-5167-02	鉗子立カバー 2号 BMC-8960	930	1,030			1186						
0-5167-03	鉗子立カバー 3号 BMC-8960	1,070	1,180			1186						
0-5167-04	鉗子立カバー 4号 BMC-8960	1,180	1,300			1186						
0-6545-01	一重四角巾 450×450mm Ad90100-544	570	630			1236						
0-6545-02	一重四角巾 600×600mm Ad90100-566	1,050	1,160			1236						
0-6545-03	一重四角巾 700×700mm Ad90100-577	1,120	1,240			1236						
0-6545-04	一重四角巾 900×900mm Ad90100-599	1,340	1,480			1236						
0-6545-05	一重四角巾 1000×1000mm Ad90100-700	1,710	1,890			1236						
0-6545-06	一重四角巾 1100×1100mm Ad90100-711	1,770	1,950			1236						
0-6545-07	一重四角巾 1300×1350mm Ad90100-733	2,550	2,850			1236						
0-6545-08	一重四角巾 1350×1600mm Ad90100-736	3,100	3,450			1236						
0-6545-09	一重四角巾 400×500mm Ad90100-545	580	640			1236						
0-6545-10	一重四角巾 450×600mm Ad90100-546	730	810			1236						
0-6545-11	一重四角巾 900×1200mm Ad90100-692	1,610	1,780			1236						
0-6545-12	一重四角巾 900×1500mm Ad90100-695	1,930	2,130			1236						
0-6545-13	一重四角巾 900×1800mm Ad90100-698	2,250	2,480			1236						
0-6545-14	一重四角巾 1100×1200mm Ad90100-712	1,980	2,180			1236						
0-6545-15	一重四角巾 1200×1200mm Ad90100-722	2,250	2,480			1236						
0-6545-16	一重四角巾 1100×1500mm Ad90100-715	2,350	2,590			1236						
0-6545-17	一重四角巾 1200×1800mm Ad90100-728	3,150	3,500			1236						
0-6545-18	一重四角巾 1600×1800mm Ad90100-766	3,850	4,250			1236						
0-6545-19	一重四角巾 1800×1800mm Ad90100-788	4,600	5,100			1236						
0-6546-01	一重四角巾(丸穴) 600×600×φ60mm Ad90301-566	1,280	1,410			1237						
0-6546-02	一重四角巾(丸穴) 900×900×φ60mm Ad90301-599	1,500	1,650			1237						
0-6546-03	一重四角巾(丸穴) 600×600×φ90mm Ad90302-566	1,280	1,410			1237						
0-6546-04	一重四角巾(丸穴) 900×900×φ90mm Ad90302-599	1,500	1,650			1237						
0-6546-06	一重四角巾(丸穴) 450×450mm Ad90301-544	710	790			1237						
0-6546-07	一重四角巾(丸穴) 1200×1200mm φ90mm Ad90301-722	2,500	2,750			1237						
0-6546-08	一重四角巾(丸穴) 1200×1200mm φ90mm Ad90302-722	2,500	2,750			1237						
0-6550-01	ゆかた型患者衣グリーン(男女兼用) S SG-1440	4,600	5,100			943						
0-6550-02	ゆかた型患者衣グリーン(男女兼用) M SG-1440 緑	4,600	5,100			943						
0-6550-03	ゆかた型患者衣グリーン(男女兼用) L SG-1440 緑	4,600	5,100			943						
0-6550-04	ゆかた型患者衣グリーン(男女兼用) LL SG-1440 緑	4,600	5,100			943						
0-6550-05	ゆかた型患者衣グリーン(男女兼用) EL SG-1440	4,600	5,100			943						
0-7279-01	電動噴霧器 フォグマスタドライブジェット トライジェット(R) 6208	115,000	125,000	2114	2114	1064						
0-7279-02	電動噴霧器 フォグマスタ マイクロジェット 高濃度微量散布方式 マイクロジェット(R) 7401	140,000	150,000	2114	2114	1064						
0-7279-03	電動噴霧器 フォグマスタ FMJR(ジュニア)	42,000	48,000	2114	2114	1064						
I-1326-01	SPカブラ TypeA 1S-A-SUS-NBR	5,390	5,930	1905	1905							
I-1326-02	SPカブラ TypeA 2S-A-SUS-NBR	5,900	6,490	1905	1905							
I-1326-03	SPカブラ TypeA 3S-A-SUS-NBR	7,550	8,320	1905	1905							
I-1326-04	SPカブラ TypeA 4S-A-SUS-NBR	16,620	18,500	1905	1905							
I-1326-05	SPカブラ TypeA 1P-A-SUS-NBR	2,220	2,460	1905	1905							
I-1326-06	SPカブラ TypeA 2P-A-SUS-NBR	2,560	2,830	1905	1905							
I-1326-07	SPカブラ TypeA 3P-A-SUS-NBR	3,200	3,530	1905	1905							
I-1326-08	SPカブラ TypeA 4P-A-SUS-NBR	6,490	7,150	1905	1905							
I-1749-01	ナスコ・ワールバック(書き込み面あり) 118mL 500枚入 B01062	22,200	22,900	2046	2046	35						
I-1749-02	ナスコ・ワールバック(書き込み面あり) 532mL 500枚入 B01065	31,600	32,600	2046	2046	35						
I-1749-03	ナスコ・ワールバック(書き込み面あり) 710mL 500枚入 B01297	40,100	41,400	2046	2046	35						
I-1749-07	ナスコ・ワールバック(書き込み面あり) 207mL 500枚入 B01489	25,400	26,200	2046	2046	35						
I-1749-08	ナスコ・ワールバック(書き込み面あり) 384mL 500枚入 B01490	35,900	37,000	2046	2046	35						
I-1749-09	ナスコ・ワールバック(書き込み面あり) 710mL 500枚入 B01196	52,100	53,700	2046	2046	35						
I-1749-11	ナスコ・ワールバック(書き込み面あり) 2041mL 500枚入 B01515	71,200	73,400	2046	2046	35						
I-1749-12	ナスコ・ワールバック(書き込み面あり) 2721mL 250枚入 B01445	61,500	63,400	2046	2046	35						
I-1749-13	ナスコ・ワールバック(書き込み面あり) 3637mL 250枚入 B01446	64,500	66,500	2046	2046	35						
I-1749-15	ナスコ・ワールバック(書き込み面あり) 5441mL 100枚入 B01447	45,300	46,700	2046	2046	35						
I-1749-16	ナスコ・ワールバック(書き込み面あり) 1627mL 500枚入 B01195	72,800	75,000	2046	2046	35						
I-2005-03	強弱精密ピンセット DURAX No.1	5,300	5,370	2592	2592							
I-2724-01	デジタル吸光度計 Checker HC(遊離塩素) HI701	14,000	15,500	579	579	1073						
I-2724-02	デジタル吸光度計 Checker HC(全塩素) HI711	14,000	15,500	579	579	1073						
I-2724-03	デジタル吸光度計 Checker HC(リン酸塩) HI713	14,500	15,000	579	579							
I-2724-04	デジタル吸光度計 Checker HC(ヨウ素) HI718	15,000	16,000	579	579							
I-2724-05	デジタル吸光度計 Checker HC(鉄) HI721	15,000	16,000	579	579							
I-2724-06	デジタル吸光度計 Checker HC(残留塩素(全塩素)ULR) HI761	15,000	16,000	579	579	1073						
I-2724-07	デジタル吸光度計 Checker HC(残留塩素(全塩素)UHR) HI771	16,000	17,500	579	579	1073						
I-2724-11	デジタル吸光度計 Checker HC(HI761用試薬) HI95761-01	7,000	7,200	579	579	1073						
I-2724-12	デジタル吸光度計 Checker HC(HI771用試薬) HI95771-01	6,900	7,000	579	579	1073						
I-3226-11	全硬度計(日常防水型)用 全濃度域用試薬 HI93735-0	35,000	38,000	587	587							
I-3226-12	全硬度計(日常防水型)用 低濃度域用試薬 HI93735-00	12,500	13,000	587	587							
I-3226-13	全硬度計(日常防水型)用 中濃度域用試薬 HI93735-01	12,500	13,000	587	587							
I-3226-14	全硬度計(日常防水型)用 高濃度域用試薬 HI93735-02	12,500	13,000	587	587							
I-3226-15	全硬度計(日常防水型)用 測定用ガラスセル(4個入り) HI731331	5,600	5,900	587	587							
I-3228-11	溶存酸素計(日常防水型)用交換電極 HI76407/4F	60,000	63,000	575	575							
I-3228-12	溶存酸素計(日常防水型)用交換電極 HI76407/10F	75,000	80,000	575	575							
I-3228-13	溶存酸素計(日常防水型)用交換電極内部液 HI7041S	5,000	5,200	575	575							
I-3250-01	ピンセット(PTFE製) 027.100 先端平 027.100(先端平)	5,200	5,360	1833	1833							
I-3250-02	ピンセット(PTFE製) 027.150 先端平 027.150(先端平)	8,300	8,560	1833	1833							
I-3538-21	pH計 本体 青 HI9811-51	50,000	52,000	564	564							
I-3538-22	pH計 本体 青 HI9812-51	51,000	53,000	564	564							
I-4315-02	低吸着チューブ 1.5mL 1袋(5000本入) T330-7LST	29,000	59,100	1107	1107							
I-4315-03	低吸着チューブ 2mL 1袋(5000本入) T330-7LST	45,000	59,900	1107	1107							
I-4458-01	イオン交換水製造装置 Milli-Diキット 卓上スタンド付 2FDJSTDKT	302,000	349,000	1461	1461							
I-5019-01	電子純水テスター PURE	19,000	21,000	573	573							
I-5019-11	電子純水テスター 用校正液 HI7033L	3,500	3,700	573	573							
I-5328-11	pHメーター pHep	10,500	10,800	568	568							
I-5329-31	ポータブルpHメーター Checker Plus用交換用電極 HI1271	5,400	5,500	568	568						506	
I-5409-01	クロスヘッド回転子ダブル 8×φ10mm 001.1110	390	370	253	253							
I-5409-02	クロスヘッド回転子ダブル 10×φ14mm 001.1114	350	410	253	253							
I-5409-05	クロスヘッド回転子ダブル 12×φ30mm 001.1130	850	880	253	253							
I-5409-06	クロスヘッド回転子ダブル 14×φ40mm 001.1140	1,300	1,340	253	253							
I-5411-01	トライアングル回転子 20×8mm 001.420	330	340	252	252							
I-5411-02	トライアングル回転子 25×8mm 001.4250	350	370	252	252							
I-5411-03	トライアングル回転子 35×9mm 001.4350	450	470	252	252							
I-5411-04	トライアングル回転子 40×14mm 001.440	560	580	252	252							
I-5411-05	トライアングル回転子 50×12mm 001.450	690	710	252	252							
I-5411-06	トライアングル回転子 55×14mm 001.4550	780	810	252	252							
I-5501-19	ベリスタポンプ 4流路 5~1500ml/h 1221334/AC-2110E-4	282,000	342,000	363	363							
I-5531-01	アイソフリーズPCRラック 紫(冷却時)・ピンク(常温時)1箱(2個入) 5640-T4	18,000	18,700	1140	1140							
I-5531-02	アイソフリーズPCRラック 緑(冷却時)・黄色(常温時)1箱(2個入) 5640-T6	18,000	18,700	1140	1140							
I-5664-21~20	酸素モニタ(OXYMAN) センサー内蔵型 OM-25MF01	100,500	105,500	616	616							
I-5692-02	DCモーター真空ポンプ(コンプレッサー兼用タイプ) DP0140-DC12V	13,930	16,030	391	391							
I-5692-03	DCモーター真空ポンプ(コンプレッサー兼用タイプ) DP0102-DC12V	23,300	26,900	391	391							
I-5692-04	DCモーター真空ポンプ(コンプレッサー兼用タイプ) DP0105-X1-DC12V	17,330	20,000	391	391							
I-5789-12	Milli-Di 最終フィルター-取付アダプター MDI02ADAP	10,100	11,700	1461	1461							
I-6510-11	日常防水型導電率計 DiST 5	25,500	26,500	571	571							
I-6510-12	日常防水型導電率計 DiST 6	25,500	26,500	571	571							
I-6511-11	日常防水pH計 pHep 4	22,500	23,500	568	568							
I-6511-12	日常防水pH計 pHep 5	25,500	26,500	568	568							
I-667-01	ハンドポンプ 本体 J7090-11	38,000	39,900	376	376</							

品番	商品名	旧価格	新価格	研究用総合 編纂2025 頁	サンクアスト 2025 頁	ナビス2026 頁	サニーフーズ 2024 頁	設備2024 頁	アズビエア 2025 頁	アズツール 2021 頁	サイエンスカ タログ 頁	災害対策カ タログ 頁
1-7397-01	ビペットチップ 200本入 バルクパック入り 4530-00	8,380	6,880	1177	1177							
1-7731-02	PTFE攪拌棒(アンカー型) φ6×400mm 005.340.2	7,100	7,320	212	212							
1-7731-04	PTFE攪拌棒(アンカー型) φ8×300mm 005.330.8.2	8,900	9,170	212	212							
1-7731-05	PTFE攪拌棒(アンカー型) φ8×400mm 005.340.8.2	9,300	9,590	212	212							
1-7731-06	PTFE攪拌棒(アンカー型) φ8×500mm 005.350.8.2	10,200	10,600	212	212							
1-7732-01	PTFE攪拌棒(プロベラ型) φ6×300mm 005.230.2	11,300	11,700	212	212							
1-7732-02	PTFE攪拌棒(プロベラ型) φ6×400mm 005.240.2	13,300	13,800	212	212							
1-7732-03	PTFE攪拌棒(プロベラ型) φ6×500mm 005.250.2	13,600	14,100	212	212							
1-7732-04	PTFE攪拌棒(プロベラ型) φ8×300mm 005.230.8.2	14,200	14,700	212	212							
1-7732-05	PTFE攪拌棒(プロベラ型) φ8×400mm 005.240.8.2	14,800	15,300	212	212							
1-7732-06	PTFE攪拌棒(プロベラ型) φ8×500mm 005.250.8.2	15,000	15,500	212	212							
1-7733-01	PTFE攪拌羽根 羽根スクエア型 52×14mm 002.052.1	770	800	213	213							
1-7733-02	PTFE攪拌羽根 羽根スクエア型 76×19mm 002.076.1	950	980	213	213							
1-7733-03	PTFE攪拌羽根 羽根スクエア型 90×28mm 002.090.1	1,200	1,240	213	213							
1-7733-04	PTFE攪拌羽根 羽根ラウンド型 65×25mm 002.1065.1	1,100	1,140	213	213							
1-7733-05	PTFE攪拌羽根 羽根ラウンド型 75×25mm 002.1075.1	1,200	1,240	213	213							
1-7733-06	PTFE攪拌羽根 羽根ラウンド型 105×25mm 002.1105.1	1,500	1,560	213	213							
1-7733-07	PTFE攪拌羽根 羽根ラウンド型 125×25mm 002.1125.1	1,500	1,550	213	213							
1-7733-08	PTFE攪拌羽根 羽根ラウンド型 150×25mm 002.1150.1	1,800	1,850	213	213							
1-7733-10	PTFE攪拌シャフト φ6×400mm 005.540.1	5,700	5,880	213	213							
1-7733-11	PTFE攪拌シャフト φ6×500mm 005.550.1	6,300	6,500	213	213							
1-7733-13	PTFE攪拌シャフト φ8×400mm 005.540.8.1	6,800	7,010	213	213							
1-7733-14	PTFE攪拌シャフト φ8×500mm 005.550.8.1	7,200	7,420	213	213							
1-7734-02	PTFE攪拌棒 φ8×100mm 003.1100	720	750	214	214							
1-7734-03	PTFE攪拌棒 φ8×150mm 003.1150	900	940	214	214							
1-7734-04	PTFE攪拌棒 φ8×200mm 003.1200	1,150	1,190	214	214							
1-7734-05	PTFE攪拌棒 φ8×250mm 003.1250	1,250	1,290	214	214							
1-7734-06	PTFE攪拌棒 φ8×300mm 003.1300	1,350	1,400	214	214							
1-7735-01	PTFE攪拌棒 鉄芯型 φ6×100mm 003.100	1,050	1,110	214	214							
1-7735-02	PTFE攪拌棒 鉄芯型 φ6×150mm 003.150	1,150	1,210	214	214							
1-7735-03	PTFE攪拌棒 鉄芯型 φ6×200mm 003.200	1,250	1,320	214	214							
1-7735-05	PTFE攪拌棒 鉄芯型 φ6×300mm 003.300	1,700	1,790	214	214							
1-7735-07	PTFE攪拌棒 鉄芯型 φ6×400mm 003.400	2,250	2,370	214	214							
1-7736-01	PTFE攪拌羽根(着脱可能型) 4枚羽根 φ6 005.606040	9,600	9,900	213	213							
1-7736-02	PTFE攪拌羽根(着脱可能型) 4枚羽根 φ8 005.608040	11,300	11,700	213	213							
1-7736-03	PTFE攪拌羽根(着脱可能型) 2枚羽根 φ6 005.708070	7,900	8,150	213	213							
1-7736-04	PTFE攪拌羽根(着脱可能型) 2枚羽根 φ8 005.708070	8,600	8,870	213	213							
1-7742-01	PTFEロータ φ33mm 039.030	2,500	2,580	1795	1795							
1-7742-02	PTFEロータ φ54mm 039.050	3,900	4,030	1795	1795							
1-7742-03	PTFEロータ φ75mm 039.070	6,900	7,110	1795	1795							
1-7742-04	PTFEロータ φ105mm 039.100	14,300	14,800	1795	1795							
1-7742-05	PTFEロータ φ156mm 039.150	25,800	26,600	1795	1795							
1-7879-01	ポータブル濁度計(キットタイプ) HI93703C	228,000	235,000	594	594							
1-7965-01	卓上型吸引システム HK-435A	105,800	121,800	385	385							
1-8097-01	キューブカブラ SPC-04PH-IVR	540	630	1901	1901							
1-8097-02	キューブカブラ SPC-06PH-IVR	540	630	1901	1901							
1-8097-03	キューブカブラ SPC-04SH-BLU	1,100	1,270	1901	1901							
1-8097-04	キューブカブラ SPC-06SH-BLU	1,100	1,270	1901	1901							
1-8547-01	スクリーキャップチューブ 0.5mL(自立型)1箱(500本入)	17,000	20,700	1114	1114							
1-8664-11	ステン柄シリコンヘラ 大 SIL-201	1,700	1,900	1841	1841							
1-8664-13	ステン柄シリコンヘラ 小 SIL-202	1,600	1,700	1841	1841							
1-9428-12	超純水製造装置Simplicity用 185nmUVランプ SYN185UV1	32,600	35,000	1460	1460							
1-9428-13	超純水製造装置Simplicity用 最終フィルター SIMFILTER	20,800	22,300	1460	1460							
1-9470-18	純水製造装置 Milli-Q用最終フィルター Biopak CDFUBIO 01	30,900	33,200	1460	1460							
1-9964-02	液量計(円盤型) 250mL 423	1,600	1,230	1739	1739							
2-1818-02	粒状活性炭(1Lポリ容器入り) 気相用4GG	5,670	7,800	1319	1319							
2-3043-71	メンブレンフィルター(セルロース混合エステル) 0.8μm×φ13mm 100枚入 AAWP01300	14,600	18,700	1241	1241							
2-3099-01	タンク残留スワブ アルルゲンスクリーニング用 1箱(60本入) ALLTEC60	18,900	19,850	1361	1361		13					
2-3775-01	滅菌希釈水 9mL 1箱(20本×20袋入) D9 PBS	17,600	18,480	1373	1373	285	39					
2-3881-31	クライオバイアル オリングシールタイプ 4mL 内ネジ 自立型 1箱(100個×10袋入) T311-4A	51,400	52,700	1115	1115							
2-3881-32	クライオバイアル オリングシールタイプ 5mL 内ネジ・丸底 滅菌済 1箱(100個×10袋入) T311-5	53,200	54,600	1115	1115							
2-3881-37	クライオバイアル用キャップインサート 青 1袋(500個入) T312-2	3,080	3,170	1115	1115							
2-3962-13	遠送管(ガスホック用)ラック HS28622	24,100	24,600	1146	1146							
2-4222-02	エアージェンレーター MAS-100 NT (TM) 1.09191.0001	1,350,000	1,390,000									
2-4226-01	サルモネラ菌測定用プレート(ベトリフィルム増地) 1箱(25枚×2袋入) 6536SALX	19,800	20,790	1356	1356							19
2-4226-02	サルモネラ菌測定用プレート(ベトリフィルム増地) 1箱(25枚×8袋入) 6537SALX	66,000	69,300	1356	1356							26
2-4226-03	サルモネラ菌測定用確認ディスク(ベトリフィルム増地) 1箱(5枚入) 6538SALX	6,600	6,930	1356	1356							26
2-4226-04	サルモネラ菌測定用確認ディスク(ベトリフィルム増地) 1箱(25枚入) 6539SALX	27,500	28,880	1356	1356							26
2-4226-11	サルモネラ菌測定用増菌基礎増地(ベトリフィルム増地) SEB500	5,200	5,460	1356	1356							26
2-4226-14	ラシボート/バリアディスタ10増地(ベトリフィルム増地) BP0288500	10,000	10,500	1356	1356							26
2-4734-03	プロビオPCRチューブ COMBI-RACK 1箱(5本入) T327-1	2,800	2,880	1098	1098							
2-4753-02	ノンワイヤラック 16mm用×30本 18748-0016	3,640	6,860	1787	1787							
2-4884-01	包埋リング 1000個入 白 M460	23,500	24,100	1435	1435	319						
2-5143-01	NMR試料管ラック Z118257	7,600	7,870	1336	1336							
2-5149-01	MEISTER ピンセット AXAL No.0-AXAL	9,400	9,500	1828	1828							
2-5149-12	MEISTER ピンセット AXAL No.4 4-AXAL	10,800	11,000	1828	1828							
2-5149-14	MEISTER ピンセット AXAL No.5A 5A-AXAL	18,100	18,300	1828	1828							
2-5149-15	MEISTER ピンセット AXAL No.5AR	11,700	11,900	1828	1828							
2-5149-17	MEISTER ピンセット AXAL No.7 7-AXAL	14,300	14,800	1828	1828							
2-5149-18	MEISTER ピンセット AXAL No.7B	14,300	14,600	1828	1828							
2-5149-41	MEISTERピンセット 115mm 7A-AXAL	12,200	12,500	1828	1828							
2-5975-01	ラビッドDPD試薬 25g 36540-33	5,300	5,600	584	584							
2-5975-03	ラビッドDPD試薬 100錠 1箱(100個入) 36541-96	4,800	5,000	584	584							
2-5975-04	ラビッドDPD試薬 100錠 1箱(100個入) 36545-96	4,800	5,000	584	584							
2-5975-12	ラビッドDPD試薬II 100包(100mg分包) 36542-97	2,400	2,500	584	584							
2-6046-01	微量透析キット(3.5kDa/10~100uL) 8連×1本 40071	7,910	8,070	1263	1263							
2-6046-02	微量透析キット(3.5kDa/10~100uL) 8連×12本 40072	94,900	96,900	1263	1263							
2-6047-01	微量透析キット(3.5kDa/~300uL) 8連×1本 40786	9,740	9,940	1263	1263							
2-6047-02	微量透析キット(3.5kDa/~300uL) 8連×12本 40787	103,000	109,000	1263	1263							
2-6048-01	微量透析キット(6~8kDa/~810~100uL) 8連×1本 40075	7,910	8,070	1263	1263							
2-6049-01	微量透析キット(6~8kDa/~300uL) 8連×1本 40788	9,740	9,940	1263	1263							
2-6049-02	微量透析キット(6~8kDa/~300uL) 8連×12本 40789	96,100	98,100	1263	1263							
2-6050-01	微量透析キット(12~14kDa/10~100uL) 8連×1本 40077	7,910	8,070	1263	1263							
2-6050-02	微量透析キット(12~14kDa/10~100uL) 8連×12本 40078	86,500	88,300	1263	1263							
2-6051-01	微量透析キット(12~14kDa/~300uL) 8連×1本 40790	9,740	9,940	1263	1263							
2-6051-02	微量透析キット(12~14kDa/~300uL) 8連×12本 40791	93,300	98,100	1263	1263							
2-6154-11	タンク残留測定スワブインサート 1箱(50本入) C1F50	8,400	8,820	1361	1361		13					
2-6371-01	フラットワイヤーサンプリングバッグ 75×185mm 1箱(500枚入) B01339	21,100	21,800	2046	2046		35					
2-6371-02	フラットワイヤーサンプリングバッグ 115×230mm 1箱(500枚入) B01341	26,400	27,200	2046	2046		35					
2-7089-05	水道水直結純水製造装置Direct-Q用 エアメントフィルター 2個入り TANKMPK03	4,300	4,610	1460	1460							
2-7089-06	超純水製造装置用最終フィルター-millipak Express20 MPGP02001	18,700	20,100	1457	1457							
2-7100-03	超純水製造装置用SynergyPak1 SYPK05A1	93,900	101,000	1461	1461							
2-743-02	カルレッツ(R)オリング P-4	5,850	9,190	1979	1979							
2-743-03	カルレッツ(R)オリング P-5	5,930	9,320	1979	1979							
2-743-04	カルレッツ(R)オリング P-6	5,970	9,380	1979	1979							
2-743-05	カルレッツ(R)オリング P-7	6,040	9,490	1979								

2026年3月1日より実施

品番	商品名	旧価格	新価格	研究用総合 編2025 頁	サンクアスト 2025 頁	ナビス2026 頁	サニーフー ズ2024 頁	設備2024 頁	アズビエア 2025 頁	アズツール 2021 頁	サイエンスカ タログ 頁	災害対策カ タログ 頁
2-743-16	カルレッツ(R)リング P-21	9,550	15,000	1979	1979							
2-743-17	カルレッツ(R)リング P-22	9,790	15,400	1979	1979							
2-7437-01	ふき取りセット ELプレート用 1箱(50セット入)6451	20,000	21,000	1368	1368		27					
2-7438-01	ベトリフィルム(TM)培地 ELプレート リステリア環境微生物測定用 50枚 6447EL	29,930	30,380	1355	1355		27					
2-7438-02	ベトリフィルム(TM)培地 ELプレート リステリア環境微生物測定用 200枚 6448EL	86,900	91,250	1355	1355		27					
2-7757-01	クリールーム用消火器 PWE-3S	39,000	42,000	2363	2363							
2-8028-01	MEISTER ピンセット SA(耐酸鋼)製 No.0 0-SA	5,300	5,480	1828	1828							
2-8028-03	MEISTER ピンセット SA(耐酸鋼)製 No.00C 00C-SA	4,900	5,060	1828	1828							
2-8028-04	MEISTER ピンセット SA(耐酸鋼)製 No.00B 00B-SA	6,700	6,930	1828	1828							
2-8028-05	MEISTER ピンセット SA(耐酸鋼)製 No.0C9 0C9-SA	4,800	4,960	1828	1828							
2-8028-06	MEISTER ピンセット SA(耐酸鋼)製 No.0C11 0C11-SA	5,000	5,170	1828	1828							
2-8028-07	MEISTER ピンセット SA(耐酸鋼)製 No.1 1-SA	4,400	4,540	1828	1828							
2-8028-08	MEISTER ピンセット SA(耐酸鋼)製 No.2 2-SA	5,000	5,170	1828	1828							
2-8028-09	MEISTER ピンセット SA(耐酸鋼)製 No.2A 2A-SA	4,700	4,870	1828	1828							
2-8028-10	MEISTER ピンセット SA(耐酸鋼)製 No.3 3-SA	5,000	5,170	1828	1828							
2-8028-11	MEISTER ピンセット SA(耐酸鋼)製 No.3C 3C-SA	5,000	5,170	1828	1828							
2-8028-12	MEISTER ピンセット SA(耐酸鋼)製 No.4 4-SA	5,200	5,360	1828	1828							
2-8028-13	MEISTER ピンセット SA(耐酸鋼)製 No.5 5-SA	6,100	6,300	1828	1828							
2-8028-14	MEISTER ピンセット SA(耐酸鋼)製 No.5A 5A-SA	6,400	6,610	1828	1828							
2-8028-15	MEISTER ピンセット SA(耐酸鋼)製 No.5AR 5AR-SA	6,600	6,820	1828	1828							
2-8028-16	MEISTER ピンセット SA(耐酸鋼)製 No.6 6-SA	6,000	6,320	1828	1828							
2-8028-17	MEISTER ピンセット SA(耐酸鋼)製 No.7 7-SA	6,800	7,010	1828	1828							
2-8028-18	MEISTER ピンセット SA(耐酸鋼)製 No.7B 7B-SA	6,500	6,700	1828	1828							
2-8028-19	MEISTER ピンセット SA(耐酸鋼)製 No.AA AA-SA	4,100	4,330	1828	1828							
2-8028-20	MEISTER ピンセット SA(耐酸鋼)製 No.F F-SA	6,400	6,740	1828	1828							
2-8028-21	MEISTER ピンセット SA(耐酸鋼)製 No.SS SS-SA	5,100	5,270	1828	1828							
2-8028-22	マイスターピンセット(バリュースセット)1袋(5本入) 0-SA-5	21,900	22,600	1828	1828							
2-8028-23	マイスターピンセット(バリュースセット)1袋(5本入) 2A-SA-5	21,100	21,800	1828	1828							
2-8028-24	マイスターピンセット(バリュースセット)1袋(5本入) 4-SA-5	22,500	23,200	1828	1828							
2-8028-25	マイスターピンセット(バリュースセット)1袋(5本入) F-SA-5	28,100	29,000	1828	1828							
2-8028-28	マイスターピンセット(バリュースセット)5本入 5AR-SA-5	27,100	28,500	1828	1828							
2-8028-52	MEISTER ピンセット 5-SA-5 5本入	26,500	27,400	1828	1828							
2-8218-03	クリーンチェア 背もたれ無し リング無し TCO-7FN	18,500	24,600	844	844							
2-8218-14	クリーンチェア TCO-9LN-Z	34,600	41,600	844	844							
2-8219-21	クラフィットマスク 息がしやすいタイプ 1箱(50枚入) EF	1,500	1,650	-	-							
2-8219-22	クラフィットマスク スタンダードタイプ 1箱(50枚入) FS	3,000	3,300	-	-							
2-8219-23	クラフィットマスク ハイグレードタイプ 1箱(30枚入) FV-R	2,400	2,640	-	-							
2-8219-24	クラフィットマスク メガネ曇り防止タイプ 1箱(50枚入) EF-K	1,650	1,820	-	-							
2-8219-25	クラフィットマスク 防臭用エコムータイプ 1箱(50枚入) EC-5	5,630	5,500	-	-							
2-8220-11	クラフィットマスク クリーンエア用 エコムーフィルター 耳かけ 1包(50枚入) EFP-P	1,500	1,650	-	-							
2-8220-12	クラフィットマスク クリーンエア用 エコムーフィルター 頭かけ(オーバーヘッド) 1包(50枚入) EFH-P	3,000	2,200	-	-							
2-8220-13	クラフィットマスク クリーンエア用 スタンダードフィルター 耳かけ 防臭活性炭 1包(50枚入) AC-P	6,500	7,150	-	-							
2-8220-14	クラフィットマスク クリーンエア用 スタンダードフィルター 頭かけ 防臭活性炭 1包(50枚入) ACH-P	7,000	7,700	-	-							
2-842-11	タンク残留測定スワブ 1箱(50本入) PRO50	14,180	14,890	1361	1361		13					
2-8606-31	土壤ダイレクトpH/℃計 交換用電極 HI2923	46,000	47,000	599	599							
2-8607-21	乳製品・半固形食品用pH/℃計 HI99161D	116,000	120,000	566	566							
2-8917-11	遊離塩素粉末試薬 1000回分 1箱(100セット入) HI93701-01	5,500	5,900	-	-	1073	91					
2-8917-12	遊離塩素粉末試薬 300回分 1箱(300セット入) HI93701-03	15,500	16,500	-	-	1073	91					
2-8917-14	全塩素粉末試薬 1000回分 1箱(100セット入) HI93711-01	5,500	5,900	-	-	1073	91					
2-8917-15	全塩素粉末試薬 300回分 1箱(300セット入) HI93711-03	15,500	16,500	-	-	1073	91					
2-8917-21	ポータブルデジタル残留塩素計 遊離塩素 HI97701C	92,000	95,000	-	-		91					
2-8917-22	ポータブルデジタル残留塩素計 遊離塩素 HI97711C	97,000	99,000	-	-		91					
2-8917-33	遊離塩素錠剤試薬 500回分 1箱(500セット入) HI93701-FJ	15,500	16,500	-	-		91					
2-8917-36	全塩素錠剤/液体試薬 300回分 1箱(300セット入) HI93701-TJ	14,000	15,000	-	-		91					
2-8917-40	ポータブルデジタル残留塩素計 遊離塩素LR/全塩素LR用/バリエーター	41,000	45,000	-	-		91					
2-8917-41	デジタル吸光度計HI 971シリーズ ガラスセル用キャップ HI 731336N 1セット(4個入)	4,000	4,300	587	587		91					
2-9140-11	秤量瓶 40×40	1,500	1,600	1708	1708							
2-9140-14	秤量瓶 50×30	1,700	1,800	1708	1708							
2-9140-15	秤量瓶 50×40	1,700	1,800	1708	1708							
2-9140-16	秤量瓶 50×50	1,700	1,800	1708	1708							
2-9140-17	秤量瓶 50×70	1,700	1,800	1708	1708							
2-9140-18	秤量瓶 50×80	1,900	2,000	1708	1708							
2-9566-02	ナイロンメッシュ 1000um	10,800	12,600	1964	1964							
2-9566-03	ナイロンメッシュ 500um	10,800	12,600	1964	1964							
2-9566-04	ナイロンメッシュ 250um	10,800	12,600	1964	1964							
2-9566-05	ナイロンメッシュ 100um	14,400	16,800	1964	1964							
2-9566-06	ナイロンメッシュ 1um	55,600	64,100	1964	1964							
2-9567-02	PETメッシュ 1000um	10,700	12,500	1964	1964							
2-9567-03	PETメッシュ 500um	10,700	12,500	1964	1964							
2-9567-04	PETメッシュ 250um	10,700	12,500	1964	1964							
2-9567-05	PETメッシュ 100um	14,300	16,700	1964	1964							
2-9567-06	PETメッシュ 1um	61,300	70,800	1964	1964							
2-9568-01	PPメッシュ 3360um	22,700	26,400	1964	1964							
2-9568-02	PPメッシュ 1000um	21,000	24,400	1964	1964							
2-9568-03	PPメッシュ 500um	11,000	12,900	1964	1964							
2-9568-04	PPメッシュ 250um	17,700	20,600	1964	1964							
2-9568-05	PPメッシュ 105um	13,700	16,000	1964	1964							
2-9569-01	フッ素樹脂メッシュ 2000um	86,400	99,600	1964	1964							
2-9569-03	フッ素樹脂メッシュ 500um	71,000	82,000	1964	1964							
2-9569-04	フッ素樹脂メッシュ 250um	55,600	64,200	1964	1964							
2-9569-06	フッ素樹脂メッシュ 70um	74,600	86,100	1964	1964							
2-9570-01	PEEKメッシュ 1000um	351,000	404,000	1964	1964							
2-9570-03	PEEKメッシュ 220um	60,400	69,700	1964	1964							
2-9570-04	PEEKメッシュ 100um	81,000	93,400	1964	1964							
2-9570-05	PEEKメッシュ 35um	75,500	87,000	1964	1964							
2-9880-13	pH・EC・DOメーター (edge) 用交換pH複合電極 HI12300	35,000	39,000	558	558							
2-9880-14	pH・EC・DOメーター (edge) 用交換pH複合電極 HI12301	38,000	42,000	558	558							
2-9880-18	pH・EC・DOメーター (edge) 用交換DO電極 HI764080	92,000	95,000	558	558							
3-1000-01	Arduino Uno アルデュイーノ A000066	6,000	7,910	2068	2068							
3-113-01	乾燥用ラック P2032	16,800	17,400	2146	2146							
3-1685-01	自立式ワールバック 75×185×0.064mm 1箱(500枚入) B01364	24,200	25,000	2046	2046							
3-1685-02	自立式ワールバック 115×230×0.076mm 1箱(500枚入) B01365	36,700	37,900	2046	2046							
3-288-03	顕微鏡スライドボックス 赤 82024-610	1,530	1,560	1452	1452							
3-293-01	ガラス器具処理ボックス 6個入 56617-804	10,080	10,400	2183	2183							
3-295-09	デイスボ丸底チューブ(虚口) 125本×4箱入 47729-586	56,300	61,500	1762	1762						585	
3-3067-01	ポリエチレンメッシュ #200 PE92-#200	3,980	4,700	1963	1963							
3-3067-02	ポリエチレンメッシュ #180 PE92-#180	3,840	4,540	1963	1963							
3-3067-03	ポリエチレンメッシュ #150 PE92-#150	3,580	4,240	1963	1963							
3-3067-04	ポリエチレンメッシュ #120 PE92-#120	2,400	2,840	1963	1963							
3-3067-05	ポリエチレンメッシュ #100 PE92-#100	1,980	2,340	1963	1963							
3-3067-06	ポリエチレンメッシュ #80 PE92-#80	1,860	2,200	1963	1963							
3-3067-07	ポリエチレンメッシュ #70 PE92-#70	1,720	2,040	1963	1963							
3-3067-08	ポリエチレンメッシュ #60 PE92-#60	1,600	1,900	1963	1963							
3-3067-09	ポリエチレンメッシュ #50 PE92-#50	1,720	2,040	1963	1963							
3-3067-10	ポリエチレンメッシュ #40 PE92-#40	1,480	1,750	1963	1963							
3-3067-11	ポリエチレンメッシュ #30 PE92-#30	1,480	1,760	1963	1963							
3-3067-12	ポリエチレンメッシュ #24 PE92-#24	1,860	2,200	1963	1963							
3-3067-13	ポリエチレンメッシュ #18 PE92-#18	1,480	1,750	1963	1963							
3-3067-14	ポリエチレンメッシュ #14 PE92-#14	2,660	3,140	1963	1963							
3-3067-15	ポリエチレンメッシュ #8 PE92-#8	2,660	3,150	1963	1963							
3-3068-01	ポリエチレンメッシュ #200 PE122-#200	4,900	5,790	1963	1963							
3-3068-02	ポリエチレンメッシュ #180 PE122-#180	4,380	5,170	1963	1963							
3-3068-03	ポリエチレンメッシュ #150 PE122-#150	4,370										

2026年3月1日より実施

品番	商品名	旧価格	新価格	研究用総合 編纂2025 頁	サンクアスト 2025 頁	ナビス2026 頁	サニーフーズ 2024 頁	設備2024 頁	アズビエア 2025 頁	アズツール 2021 頁	サイエンスカ タログ 頁	災害対策カ タログ 頁
4-4845-02	カンホーグB5基礎培地 G398-10L	3,200	3,460	1414	1414							
4-4846-01	ムンダグスクグ培地 M519-1L	1,300	1,420	1413	1413							
4-4846-02	ムンダグスクグ培地 M519-10L	3,500	3,790	1413	1413							
4-4848-01	植物培養容器 C590-12EA 12個入	11,500	12,500	1097	1097							
4-4848-02	植物培養容器 C1892-12EA 12個入	13,100	13,400	1097	1097							
4-4848-04	植物培養容器 C170-10EA 10個入	13,100	14,200	1097	1097							
4-4848-05	植物培養容器 C170-25EA 25個入	21,300	23,000	1097	1097							
4-4998-02	スピナーアスコセルスター 356882	127,000	130,000	1094	1094							
4-4998-03	スピナーアスコセルスター 356884	133,000	145,000	1094	1094							
4-5104-11	クリーンルーム用カーテータブ 1160xx050-P3D-WH	3,200	3,770	2509	2509							
4-5104-12	クリーンルーム用カーテータブ 1160xx050-P3D-RE	3,200	3,770	2509	2509							
4-5104-13	クリーンルーム用カーテータブ 1160xx050-P3D-YE	3,200	3,770	2509	2509							
4-5104-15	クリーンルーム用カーテータブ 1160xx100-P3D-WH	4,200	4,940	2509	2509							
4-5104-16	クリーンルーム用カーテータブ 1160xx100-P3D-RE	4,200	4,940	2509	2509							
4-5104-17	クリーンルーム用カーテータブ 1160xx100-P3D-YE	4,200	4,940	2509	2509							
4-5436-01	メッシュカーボン510μ	21,100	24,400	1964	1964							
4-5436-02	メッシュカーボン245μ	21,100	24,400	1964	1964							
4-5436-03	メッシュカーボン105μ	26,600	30,800	1964	1964							
4-5485-02	クリーンルームスワブ LTO70PESD20	67,600	72,200						157			
4-5722-05	アクリルパイプ φ40×5×1000mm PMMA-Φ40-5	8,270	8,680	1960	1960							
4-5722-07	アクリルパイプ φ50×5×1000mm PMMA-Φ50-5	9,480	10,100	1960	1960							
4-6073-01	手動高耐荷重XYステージ(ダイヤルロック付) SHXY150-70-L	250,000	160,000	807	807							
4-6073-02	手動高耐荷重XYステージ(ダイヤルロック付) SHXY200-70-L	292,000	179,000	807	807							
4-6073-03	手動高耐荷重XYステージ(ダイヤルロック付) SHXY250-60-L	305,000	199,000	807	807							
4-6074-01	手動高耐荷重XYステージ(ダイヤルロック付) SHXY150-70-DL	211,000	184,000	807	807							
4-6074-02	手動高耐荷重XYステージ(ダイヤルロック付) SHXY200-70-DL	245,000	208,000	807	807							
4-6074-03	手動高耐荷重XYステージ(ダイヤルロック付) SHXY250-60-DL	256,000	216,000	807	807							
4-6268-04	エロジーLDPEボトル 0303-0100	440	480	1518	1518							
4-6547-01	バイオハザードゴミ箱 BowTie 本体 30L A8000B	21,800	22,400	2166	2166							
4-6618-05	ビーカー KJMAX(R) 100mL Φ50×70mm 14000-100	940	970	1716	1716							
4-6624-04	ネジロメディウム瓶 GL45 1,000mL 14395-1000	3,430	3,530	1584	1584							
4-6724-04	PTFEジョイントスリーブ 020.124	2,700	2,790	1686	1686							
4-6724-05	PTFEジョイントスリーブ 020.129	3,500	3,610	1686	1686							
4-6728-01	固定タプチューブフローター PREPSAFE アタッチメント付 12本掛け ブルー F18742-0125	10,800	11,700	1141	1141							
5-3500-21	タップ式ふるい振とう機 防音カバー無 AS200tap(50Hz)	920,000	1,050,000	323	323							
5-3500-22	タップ式ふるい振とう機 防音カバー無 AS200tap(60Hz)	920,000	1,050,000	323	323							
5-3500-23	タップ式ふるい振とう機 防音カバー付 AS200tap(50Hz)	1,420,000	1,630,000	323	323							
5-3500-24	タップ式ふるい振とう機 防音カバー付 AS200tap(60Hz)	1,420,000	1,630,000	323	323							
5-4001-13	電磁式ふるい振とう機 AS200 DigiCo(50/60Hz)	700,000	770,000	323	323							
5-4001-17	電磁式ふるい振とう機 AS200 Control(50/60Hz)	850,000	940,000	323	323							
5-5043-01	フッ素樹脂成型ひしゃく 036.100	15,500	16,000	1847	1847							
5-5043-02	フッ素樹脂成型ひしゃく 036.250	18,900	19,500	1847	1847							
5-5043-03	フッ素樹脂成型ひしゃく 036.500	24,100	24,900	1847	1847							
5-5043-04	フッ素樹脂成型ひしゃく 036.1000	32,200	33,200	1847	1847							
5-5043-11	フッ素樹脂成型ひしゃく用 延長用ねじ込取手 036.001	10,500	10,900	1847	1847							
5-5366-11	フッ素樹脂製三方コック(ロック締め) 016.1302.6.2	6,000	6,180	1894	1894							
5-5367-12	フッ素樹脂製三方コック(ロック締め) 016.1204.2	4,800	4,950	1894	1894							
5-5368-11	フッ素樹脂製三方コック(差し込み止) 016.802.5.2	4,200	4,330	1894	1894							
5-5368-12	フッ素樹脂製三方コック(差し込み止) 016.802.2	4,600	4,740	1894	1894							
5-5368-13	フッ素樹脂製三方コック(差し込み止) 016.803.2	5,000	5,150	1894	1894							
5-5368-14	フッ素樹脂製三方コック(差し込み止) 016.804.2	5,500	5,670	1894	1894							
5-5369-11	フッ素樹脂製三方コック(差し込み止) 016.7025.2 016.702.5.2	3,900	4,020	1894	1894							
5-5369-12	フッ素樹脂製三方コック(差し込み止) 016.702.2	3,900	4,030	1894	1894							
5-5369-13	フッ素樹脂製三方コック(差し込み止) 016.703.2	4,100	4,230	1894	1894							
5-5369-14	フッ素樹脂製三方コック(差し込み止) 016.704.2	4,100	4,230	1894	1894							
5-5600-11	電磁ふるい振とう機 湿式用頭部蓋(ふるい用) 31.0400.00	190,000	195,000	322	322							
5-5600-21	ふるい振とう機乾式用受皿(ふるい用)1	24,000	25,000	322	322							
5-5600-41	電磁ふるい振とう機 AC100~120V spartan	574,000	641,000	322	322							
5-5600-42	電磁ふるい振とう機 AC100~120V pro	702,000	796,000	322	322							
5-5600-43	電磁ふるい振とう機 AC100~120V BASIC	553,000	620,000	322	322							
5-5600-51	電磁ふるい振とう機 AC200~240V spartan	574,000	641,000	322	322							
5-5600-52	電磁ふるい振とう機 AC200~240V pro	702,000	796,000	322	322							
5-5600-53	電磁ふるい振とう機 AC200~240V BASIC	553,000	620,000	322	322							
5-5663-11	ステンレスコーキングヘラ 10×97×190	690	900	1842	1842							
5-5663-12	ステンレスコーキングヘラ 15×115×210	860	1,150	1842	1842							
5-5663-13	ステンレスコーキングヘラ 24×145×250	990	1,350	1842	1842							
60-9607-71	導電性両面粘着テープ 200×300mm No.118M	2,840	1,070	-	-							
61-0120-24	デジタル糖度計 HI96801	45,000	49,000	704	704		102					
61-0120-25	デジタル屈折計 HI96802	48,000	58,000	704	704		102					
61-0120-26	デジタル屈折計 HI96803	48,000	58,000	704	704		102					
61-0120-27	デジタル屈折計 HI96804	48,000	58,000	704	704		102					
61-0120-28	デジタル屈折計 HI96812	48,000	58,000	704	704		102					
61-0120-30	デジタル糖度計 HI96813	48,000	58,000	704	704		102					
61-0120-31	デジタル屈折計 HI96814	48,000	58,000	704	704		102					
61-0120-32	デジタル塩分計(食塩用) HI96821	45,000	49,000	704	704		102					
61-0120-33	デジタル塩分計(海水用) HI96822	45,000	49,000	704	704							
61-0120-34	デジタル屈折計 HI96831	48,000	58,000	704	704							
61-0120-35	デジタル屈折計 HI96832	48,000	58,000	704	704							
61-0120-62	チェックテンブ 温度デスター HI98501	10,400	11,000	512	512							
61-0120-63	チェックテンブ1 温度デスター HI98509	10,800	12,000	512	512							
61-0757-81	術前術後衣 S グリーン EG315	8,600	9,500			944						
61-0757-82	術前術後衣 M グリーン EG315	8,600	9,500			944						
61-0757-83	術前術後衣 L グリーン EG315	8,600	9,500			944						
61-3640-06	環境対応OATapp 3P 4個口 3m T-ECOH3430NM	2,760	3,780									
61-4385-64	2方ニップル電磁弁 FSS-03TBFC DQ24V F-2361-003	11,300	14,100	1895	1895							
61-4385-65	2方ニップル電磁弁 FSS-03TBFC AC100V F-2361-004	14,500	17,000	1895	1895							
62-0872-89	セラチューブ(ブリストルバック) イエロー #1TB-11	1,600	1,700			680						
62-9201-36	ビタミンシュース(メンズ) ホワイト 24 MN415	4,800	5,300			1956						
62-9201-37	ビタミンシュース(メンズ) ホワイト 24.5 MN415	4,800	5,300			1956						
62-9201-38	ビタミンシュース(メンズ) ホワイト 25 MN415	4,800	5,300			1956						
62-9201-39	ビタミンシュース(メンズ) ホワイト 25.5 MN415	4,800	5,300			1956						
62-9201-40	ビタミンシュース(メンズ) ホワイト 26 MN415	4,800	5,300			1956						
62-9201-41	ビタミンシュース(メンズ) ホワイト 26.5 MN415	4,800	5,300			1956						
62-9201-42	ビタミンシュース(メンズ) ホワイト 27 MN415	4,800	5,300			1956						
62-9201-43	ビタミンシュース(メンズ) ホワイト 28 MN415	4,800	5,300			1956						
63-1269-48	手動高耐荷重XYステージ SHXY150-150	170,000	140,000	807	807							
63-1269-49	手動高耐荷重XYステージ SHXY200-200	200,000	170,000	807	807							
63-1269-50	手動高耐荷重XYステージ SHXY250-200	221,000	185,000	807	807							
63-1269-58	薄型・精密・大型XYステージ(マイクロヘッド駆動) 高耐荷重タイプ IXY125	260,000	150,000	807	807							
63-1269-59	薄型・精密・大型XYステージ(マイクロヘッド駆動) 高耐荷重タイプ IXY150	313,000	220,000	807	807							
63-1269-60	薄型・精密・大型XYステージ(マイクロヘッド駆動) 高耐荷重タイプ IXY175	337,000	230,000	807	807							
63-1269-61	薄型・精密・大型XYステージ(マイクロヘッド駆動) 高耐荷重タイプ IXY200	364,000	250,000	807	807							
63-1536-79	GroLine 土壌ダイレクタpHセンサー(Soil pH Tester) H981030	26,000	29,000	600	600							
63-1536-84	食肉用pH/℃計 HI99163D	116,000	120,000				89					
63-1829-35	ふき取り用ドライスポンジ 100個入 BP237SPG	35,000	38,500				13					
63-1829-36	ふき取り用ドライスポンジスティック 100個入 SSL100	37,000	40,700				13					
63-1829-37	ふき取り用中和緩衝液含浸スポンジ 100個入	36,000	50,400				13					
63-1829-38	ふき取り用中和緩衝液含浸スポンジスティック 100個入 SSL10NB	38,000	53,200				13					
63-2744-82	GroLine ホータブpH/EC/TDS/℃計 HI9814D	56,000	58,000	600	600							
63-2744-96	ワイン用pH/℃計 HI99111D	108,000	115,000				89					
63-2744-99	ビール用pH/℃計 HI99151D	108,000	115,000				89	</				

品番	商品名	旧価格	新価格	研究用総合 編纂2025 頁	サンクアスト 2025 頁	ナビス2026 頁	サニーフーズ 2024 頁	設備2024 頁	アズビエア 2025 頁	アズツール 2021 頁	サイエンスカ タログ 頁	災害対策カ タログ 頁
64-5225-12	調理用デジタル温度テスター(ホワイト:乳製品) HI151	9,500	9,800				76					
64-5225-13	調理用デジタル温度テスター(レッド:生肉) HI151-1	9,500	9,800				76					
64-5225-14	調理用デジタル温度テスター(ブルー:生魚) HI151-2	9,500	9,800				76					
64-5225-15	調理用デジタル温度テスター(イエロー:調理肉) HI151-3	9,500	9,800				76					
64-5225-16	調理用デジタル温度テスター(グリーン:サラダ/フルーツ) HI151-4	9,500	9,800				76					
64-5225-17	調理用デジタル温度テスター(ブラウン:野菜) HI151-5	9,500	9,800				76					
64-8893-10	顕厚計 TG-01U	25,500	19,500	715	715							
64-8893-11	顕厚計 TG-02U	34,600	27,000	715	715							
64-9335-07	外付けSSD ノック式 USB3.2 (Gen2)対応 250GB ブラック ESD-EPK0250GBK	9,240	12,960	2613	2613							
64-9335-10	外付けSSD ノック式 USB3.2 (Gen2)対応 500GB ブラック ESD-EPK0500GBK	14,300	19,220	2613	2613							
64-9335-13	外付けSSD ノック式 USB3.2 (Gen2)対応 1TB ブラック ESD-EPK1000GBK	23,000	32,210	2613	2613							
64-9634-87	ポータブル項目別 吸光度計(全硬度) HI97735	72,000	76,000	587	587							
65-0535-86	HALO2 ワイヤレスpHテスター 土壌 HI 9810302	30,000	32,000	568	568		89					
65-0535-87	HALO2 ワイヤレスpHテスター ビール HI 9810312	31,000	32,000	568	568		89					
65-0535-88	HALO2 ワイヤレスpHテスター チーズ HI 9810322	31,000	32,000	568	568		89					
65-0535-89	HALO2 ワイヤレスpHテスター ワイン HI 9810332	31,000	32,000	568	568		89					
65-0535-90	HALO2 ワイヤレスpHテスター ミルク HI 9810342	31,000	32,000	568	568		89					
65-0535-91	HALO2 ワイヤレスpHテスター 寿司 HI 9810352	31,000	32,000	568	568		89					
65-0535-92	HALO2 ワイヤレスpHテスター 食肉 HI 9810362	31,000	32,000	568	568		89					
65-0535-93	HALO2 ワイヤレスpHテスター 皮膚/頭皮 HI 9810372	31,000	32,000	568	568		89					
65-0535-94	HALO2 ワイヤレスpHテスター パン生地 HI 9810382	31,000	32,000	568	568		89					
65-0535-95	HALO2 ワイヤレスpHテスター チョコレート HI 9810392	31,000	32,000	568	568		89					
65-0535-96	HALO2 ワイヤレスpHテスター ラボ(KCL) HI 9810402	32,000	35,000	568	568							
65-0535-97	HALO2 ワイヤレスpHテスター ラボ(ゲル) HI 9810412	34,000	35,000	568	568							
65-0535-98	高精度pHテスター HALO2(ヘイロー2) フィールド用(プラスチックボディ) HI 9810422	31,000	32,000	568	568							
65-0535-99	HALO2 ワイヤレスpHテスター 化粧クリーム HI 9810432	34,000	35,000	568	568							
65-0536-01	HALO2 ワイヤレスpHテスター 紙/皮革 HI 9810442	34,000	35,000	568	568							
65-3411-68	リンクステーション SOHO向け LS710DNB ネットワークHDD 1ベイ 1TB LS710DN0101B	38,600	39,900	2614	2614							
65-3418-22	リンクステーション SOHO向け LS710DNB ネットワークHDD 1ベイ 2TB LS710DN0201B	42,000	42,400	2614	2614							
65-8735-71	イス付きナースカート MC-NA-CNN	281,600	361,000			1400						
65-9616-40	ゲル撮影装置 Limited-STAGE II-VP 感熱プリンター、8インチモニター付き	590,000	630,000	1302	1302							
65-9616-41	ゲル撮影装置 Limited-STAGE II-UVVP UV付き、感熱プリンター、8インチモニター付き	740,000	780,000	1302	1302							
65-9616-42	ゲル撮影装置 Limited-STAGE II-100VP 青色LED付き、感熱プリンター、8インチモニター付き	730,000	770,000	1302	1302							
66-0023-87	チップアイスマーカ 500×450×800mm CM-60A	678,000	932,000	185	185							
66-0023-88	チップアイスマーカ 600×600×800mm CM-100K	872,000	1,145,000	185	185							
66-0023-89	チップアイスマーカ 900×600×800mm CM-100K-50	948,000	1,244,000	185	185							
66-0023-90	チップアイスマーカ 1060×600×800mm CM-200K	1,244,000	1,634,000	185	185							
66-0023-91	フレックアイスマーカ 600×600×800mm FM-120K	794,000	1,043,000	186	186							
66-0023-92	フレックアイスマーカ 900×600×800mm FM-120K-50	872,000	1,145,000	186	186							
67-2479-63	ワンタッチイベントテント 切妻アルミ3060 ホワイト 6179-01	283,000	293,000			903						
67-2479-64	ワンタッチイベントテント 切妻アルミ3060 レッド 6179-10	283,000	293,000			903						
67-2479-65	ワンタッチイベントテント 切妻アルミ3060 イエロー 6179-30	283,000	293,000			903						
67-2479-66	ワンタッチイベントテント 切妻アルミ3060 グリーン 6179-40	283,000	293,000			903						
67-2479-67	ワンタッチイベントテント 切妻アルミ3060 ブルー 6179-51	283,000	293,000			903						
67-2479-68	ワンタッチイベントテント 切妻アルミ3060 ブラック 6179-99	283,000	293,000			903						
67-4597-20	電位差測定装置 HI 931-01	882,000	934,000	574	574							
67-5757-86	携帯型暑さ/乾燥指数計 HV-610	2,600	1,700	658	658		167				535	
67-6640-53	HI 9811-51/ HI 9812-51用交換電極 ケーブル1m付 HI 1285-51	34,000	35,000	564	564							
67-7741-01	リニア真空ポンプ 7L/min VP0125	27,600	31,900	402	402							
67-7741-02	リニア真空ポンプ 60L/min VP0940	69,800	80,300	402	402							
67-7913-01	精密作業用はさみ 1C200/EL	8,500	8,600	2061	2061							
68-8331-01	フッ素樹脂成型センサー 500.300.1	22,200	23,400	505	505							
68-8331-02	フッ素樹脂成型センサー 600.301.1	17,800	18,700	505	505							
68-8331-03	フッ素樹脂成型センサー 500.100.1	17,400	18,300	505	505							
68-8331-04	フッ素樹脂成型センサー 500.200.1	22,300	23,500	505	505							
68-8331-05	フッ素樹脂成型センサー 600.101.1	12,500	13,200	505	505							
68-8331-06	フッ素樹脂成型センサー 600.201.1	14,400	15,200	505	505							
68-8641-16	拭き取り検査用 クイックワブ 25本×2袋 6432	9,800	10,780	1369	1369		15					
68-8641-22	ベトリフィルム(TM)培地(水中一般生菌測定用/50枚×2袋) 6450AQHC	16,300	17,210	1355	1355		26					
68-8641-23	ベトリフィルム(TM)培地(水中一般生菌測定用/50枚×20袋) 6452AQHC	115,500	121,280	1355	1355		26					
68-8746-32	護身用放電付きLEDライト a29959	4,530	3,900									276
68-8689-01	白金線 φ0.3mm×1m 1本入	65,000	80,400	1924	1924		54					
68-8689-02	白金線 φ0.5mm×1m 1本入	120,000	153,000	1924	1924		54					
68-8690-01	白金耳 φ0.3×50mm×φ2.5本入	25,900	30,400	1924	1924	314	55					
68-8690-02	白金耳 φ0.3×100mm×φ2.5本入	35,600	43,200	1924	1924	314	55					526
68-8690-03	白金耳 φ0.5×50mm×φ5.5本入	57,400	68,000	1924	1924	314	55					526
68-8690-04	白金耳 φ0.5×100mm×φ5.5本入	79,400	95,800	1924	1924	314	55					526
68-8754-03	マグネチックフィルターファンネル 300mL(フタ無し) 4242	37,900	39,000	1255	1255							
68-8762-01	ナスコ・ワールバック(書き込み面なし) 118mL 500枚入 B00679	17,900	18,500	2046	2046		35					
68-8762-02	ナスコ・ワールバック(書き込み面なし) 532mL 500枚入 B00736	25,900	26,700	2046	2046		35					
68-8762-08	ナスコ・ワールバック(書き込み面なし) 710mL 500枚入 B01020	32,700	33,700	2046	2046		35					
68-8762-13	ナスコ・ワールバック(書き込み面なし) 1627mL 250枚入 B01532	60,300	62,200	2046	2046		35					
69-3870-10	サーモキャプチャーIR THG-02	58,000	49,800									
7-149-02	MDーハンド吸着パッド φ4mm 10コ入 吸着パッドφ4mm(10個入)	3,690	4,250	2579	2579							
7-149-03	MDーハンド吸着パッド φ6mm 10コ入 吸着パッドφ6mm(10個入)	3,690	4,250	2579	2579							
7-149-04	MDーハンドピンセット P-100	11,790	13,560	2579	2579							
7-149-11	MDーハンド MH-100A	63,500	73,200	2579	2579							
7-2586-01	オレンジア 救急絆創膏 M 1箱(100枚入)	200	250			1987						
7-3860-01	ステップピン エコフライ シルバー 20L EK9377MT-20L	16,000	18,000			1097						
7-3860-02	ステップピン エコフライ シルバー 30L EK9377MT-30L	21,000	22,000			1097						
7-3860-03	ステップピン エコフライ シルバー 45L EK9377MT-45L	25,000	27,000			1097						
7-3860-04	ステップピン エコフライ シルバー 20×2L EK9377MT-20L+20L	26,000	28,000			1097						
7-4346-01	入れ歯ケース(書き込みスペース付き) ピンク H980-91	150	160			1716						
7-4346-02	入れ歯ケース(書き込みスペース付き) グリーン H980-92	150	160			1716						
7-4346-51	入れ歯ケース(書き込みスペース付き) ピンク 200個 H980-91	27,000	27,800			1716						
7-6518-01	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ50mm 2個入 700/300/050JP	2,800	2,900			709						
7-6518-02	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ53mm 2個入 700/300/053JP	2,800	2,900			709						
7-6518-03	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ56mm 2個入 700/300/056JP	2,800	2,900			709						
7-6518-04	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ59mm 2個入 700/300/059JP	2,800	2,900			709						
7-6518-05	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ62mm 2個入 700/300/062JP	2,800	2,900			709						
7-6518-06	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ65mm 2個入 700/300/065JP	2,800	2,900			709						
7-6518-07	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ68mm 2個入 700/300/068JP	2,800	2,900			709						
7-6518-08	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ71mm 2個入 700/300/071JP	2,800	2,900			709						
7-6518-09	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ74mm 2個入 700/300/074JP	2,800	2,900			709						
7-6518-10	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ77mm 2個入 700/300/077JP	2,800	2,900			709						
7-6518-11	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ80mm 2個入 700/300/080JP	2,800	2,900			709						
7-6518-12	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ85mm 2個入 700/300/085JP	2,800	2,900			709						
7-6518-13	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ90mm 2個入 700/300/090JP	2,800	2,900			709						
7-6518-14	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ95mm 2個入 700/300/095JP	2,800	2,900			709						
7-6518-15	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ100mm 2個入 700/300/100JP	2,800	2,900			709						
7-6518-16	Wallaceリングベッサリー(子宮下垂・子宮脱治療用リング) φ110mm 2個入 700/300/110JP	2,800	2,900			709						
7-652-01	パデューベグボード 300×455×20mm 12-3000	61,800	69,000			688						
7-7235-01	Xキューブステップピン 20L EK9368MT	18,000	21,000			1097						
7-7235-02	Xキューブステップピン 30L EK9368MT	23,000	25,000			1097						
7-7235-03	Xキューブステップピン 45L EK9368MT	26,000	28,000			1097						
7-7235-04	Xキューブステップピン 20+20L EK9368MT	27,000	30,000			1097						
7-773-01	哺乳びん殺菌乾燥保管											

品番	商品名	旧価格	新価格	研究用総合 編纂2025 頁	サンクアスト 2025 頁	ナビス2026 頁	サニーフー 2024 頁	設備2024 頁	アズビエア 2025 頁	アズツール 2021 頁	サイエンスカ タログ 頁	災害対策カ タログ 頁
8-3447-02	歯垢染め出し液 (FRUITYライト) グレーブ風味 5mL		430	500			1708					
8-465-02	遮光チャック袋 MGUV3P0305		830	900			1567					
8-465-04	遮光チャック袋 MGUV3P0406		960	1,050			1567					
8-465-08	遮光チャック袋 MGUV3P0814		4,740	5,150			1567					
8-4770-11	小型ニッパ式爪切り 直刃 全長100mm NAIL-SH		3,500	4,500			639					
8-4770-12	小型ニッパ式爪切り 直刃 全長130mm NAIL-LH		3,750	4,850			639					
8-4770-13	小型ニッパ式爪切り アングル刃 全長100mm NAIL-SAT		3,500	4,500			639					
8-4770-14	小型ニッパ式爪切り アングル刃 全長130mm NAIL-LAT		3,750	4,850			639					
8-494-02	プラスチックチューブ用蓋 H21011-0000		2,590	4,370	2060	2060						
8-7671-11	ゼルボックスデジタル XL(電子補償器) 10000250		76,000	81,900			1929					
89-0126-10	ナフタリン 01422-1G 01422		7,820	8,290	1500	1500						
89-0126-49	アクリルアミド溶液 01697-500ML 01697		59,500	61,500	1500	1500						
89-0127-40	塩酸 2104-50ML 2104		11,300	11,700	1500	1500						
89-0130-78	Ehrlichs solution 03891-250ML 03891		13,100	13,500	1500	1500						
89-0131-09	Acesulfame K, for food analysis, >=99.0% 04054-100G 04054		43,300	44,800	1500	1500						
89-0131-10	Acesulfame K, for food analysis, >=99.0% 04054-25G 04054		15,700	16,200	1500	1500						
89-0173-02	クロロホルム:イソアミルアルコール 24:1 25666-100ML 25666		20,200	20,900	1500	1500						
89-0173-03	クロロホルム:イソアミルアルコール 24:1 25666-500ML 25666		83,100	86,000	1500	1500						
89-0185-64	Diethylenetriaminepentaacetic acid, for complexometry, >=99.0% 32319-100G-F 32319		13,400	14,000	1500	1500						
89-0185-65	Diethylenetriaminepentaacetic acid, for complexometry, >=99.0% 32319-500G-F 32319		54,000	56,600	1500	1500						
89-0194-45	ジクロロメタン(高速液体クロマトグラフィー用) 34856-100ML 34856		4,410	4,570	1500	1500						
89-0195-21	蒸留水 34877-1L 34877		10,100	10,400	1500	1500						
89-0205-03	Melting point standard 235-237 C, analytical standard 41019-1G 41019		3,320	3,520	1500	1500						
89-0205-04	Melting point standard 235-237 C, analytical standard 41019-5G 41019		7,620	8,080	1500	1500						
89-0206-87	Melting point standard 182-184 C, analytical standard 42183-1G 42183		10,180	10,800	1500	1500						
89-0206-88	Melting point standard 182-184 C, analytical standard 42183-250MG 42183		3,030	3,220	1500	1500						
89-0210-92	Iron(III) citrate tribasic monohydrate, 18-20% Fe basis (T) 44941-250G 44941		13,000	13,600	1500	1500						
89-0210-94	Iron(II) gluconate hydrate 44948-250G 44948		14,800	15,500	1500	1500						
89-0227-19	Melting point standard 47-49 C, analytical standard 50296-1G 50296		7,540	4,490	1500	1500						
89-0227-20	Melting point standard 47-49 C, analytical standard 50296-5G 50296		18,000	17,500	1500	1500						
89-0246-17	Potassium chloride solution, BioUltra, ~3 M in H2O 60135-250ML 60135		15,400	15,900	1500	1500						
89-0246-29	Potassium citrate monobasic, >=98.0% (NT) 60215-1KG 60215		42,400	43,900	1500	1500						
89-0246-42	ケイフ化カリウム 60322-1KG 60322		16,000	16,700	1500	1500						
89-0255-65	Melting point standard 283-286 C, analytical standard 67372-1G 67372		33,500	35,600	1500	1500						
89-0255-66	Melting point standard 283-286 C, analytical standard 67372-250MG 67372		11,300	12,000	1500	1500						
89-0272-20	Melting point standard 121-123 C, analytical standard 76170-5G 76170		8,610	9,130	1500	1500						
89-0273-56	Phosphate buffer solution, BioUltra, >0.044 M Na2HPO4, >0.022 M NaH2PO4 76847-100ML-F 76847		37,200	38,500	1500	1500						
89-0274-97	フェノール:クロロホルム:イソアミルアルコール 125:24:1 77619-100ML 77619		19,300	19,900	1500	1500						
89-0274-98	フェノール:クロロホルム:イソアミルアルコール 125:24:1 77619-500ML 77619		73,500	76,000	1500	1500						
89-0275-04	Poly(vinylpyrrolidone), ~110 mum particle size 77627-25G 77627		11,000	11,400	1500	1500						
89-0277-34	Amino Acids Mix Solution, certified reference material, TraceCERT(R) 79248-5X2ML 79248		69,400	73,400	1500	1500						
89-0286-17	Spermine dihydrate, BioUltra, >=99.5% (GC) 85588-5G 85588		44,200	45,700	1500	1500						
89-0286-64	Methoxyethylene glycol 5.000 propionic acid N-succinimidyl ester, >=80% 85969-1G 85969		154,000	159,000	1500	1500						
89-0286-67	Methoxyethylene glycol 5.000 acetic acid N-succinimidyl ester, >=80% 85973-1G 85973		102,000	105,000	1500	1500						
89-0289-30	Tetrabutylammonium methanesulfonate, >=97.0% (T) 86877-10G 86877		25,600	26,400	1500	1500						
89-0301-71	TE緩衝液 pH8.0 分子生物学用 93283-100ML 93283		6,120	6,410	1490	1490						
89-0301-87	トリス(ヒドロキシメチル)アミノメタン塩酸塩、1M 標準原液 pH7.40 分子生物学用 93313-1L 93313		24,000	24,800	1500	1500						
89-0302-51	トリトンX-100溶液 93443-100ML 93443		32,300	33,900	1490	1490						
89-0307-15	Xanthinol, for the detection of urea, >=99.0% (HPLC) 95580-100G 95580		133,000	139,000	1500	1500						
89-0307-16	Xanthinol, for the detection of urea, >=99.0% (HPLC) 95580-25G 95580		40,000	41,900	1500	1500						
89-0416-75	蒸留水 270733-1L 270733		9,090	9,420	1500	1500						
89-0435-11	トリフルオロ酢酸 302031-100ML 302031		42,600	44,100	1500	1500						
89-0468-48	クロロホルム 366927-100ML 366927		12,400	12,800	1500	1500						
89-0615-59	ジクロロメタン, Methylene chloride 650463-1L 650463		14,900	15,400	1500	1500						
89-0615-63	クロロホルム 650471-1L 650471		21,700	22,400	1500	1500						
89-0615-67	クロロホルム 650498-1L 650498		22,000	22,700	1500	1500						
89-0769-81	ORBO(TM) 65P (OVS) Amberlite(R) XAD(R)-4 specially cleaned, 160/80 mg, RC.GFF.F.F (specially 20029-U		63,400	67,000	1315	1315						
89-0770-15	ORBO(TM) 32 Small Activated Coconut Charcoal (20/40), 100/50 mg, W.F.F separators, O.D. x L 20267-U		17,900	18,900	1315	1315						
89-0770-47	ORBO(TM) 1000 Precleaned Small PUF Plug, Replacement, O.D. x L 22 mm x 7.6 cm, pkg of 3 ea 20600-U		10,300	10,900	1315	1315						
89-0771-91	ORBO(TM) 1500 Precleaned Small PUF/Amberlite(R) XAD(R)-2/PUF Cartridge, O.D. x L 22 mm x 10 21233-U		35,000	37,000	1315	1315						
89-0788-87	Tenax(R) TA, 35/60 mesh, glass TD tube, fritted, O.D. x L 1/4 in. x 3 1/2 in., preconditioned 29530-U		197,000	208,000	1315	1315						
89-0788-92	Carbopak(TM) B, glass TD tube, fritted, O.D. x L 1/4 in. x 3 1/2 in., preconditioned, pkg o 29535-U		147,400	156,000	1315	1315						
89-0789-58	Tenax(R) GR, 60/80 mesh, O.D. x L 1/4 in. x 3 1/2, stainless steel TD tube, sealed with Bras 29742-U		171,000	180,000	1315	1315						
89-0789-66	Carbotrap(R) 217, O.D. x L 1/4 in. x 3 1/2 in., glass TD tube, fritted, Sealed with Brass En 29750-U		170,000	179,000	1315	1315						
89-0790-21	Tenax(R) TA, glass TD tube, fritted, O.D. x L 1/4 in. x 3 1/2 in., unconditioned, pkg of 10 30134-U		107,000	113,000	1315	1315						
89-0801-14	ORBO(TM)-5020 Tenax(R) TA (35/60) OVS Tube, pkg of 10 ea 52681-U		53,500	56,500	1315	1315						
89-0804-16	LpDNPH Rezorin(TM) Cartridge, bed wt. 350 mg, volume 3 mL, pkg of 10 ea 54074-U		19,100	20,200	1315	1315						
89-0804-17	LpDNPH Rezorin(TM) Cartridge, bed wt. 350 mg, volume 3 mL, pkg of 50 ea 54075-U		69,800	73,700	1315	1315						
89-0804-84	Carboxen(R) 572 for Tobacco Smoke, Adsorbent Bed Wt: 300 mg, cartridge size 3 mL, pkg of 50 54294-U		110,000	116,000	1315	1315						
89-0810-59	SPME fiber assembly Polydimethylsiloxane (PDMS), df 30 mum (nonbonded), needle size 23 ga, f 57289-U		63,300	67,300	1315	1315						
89-0810-60	SPME fiber assembly Polydimethylsiloxane (PDMS), df 7 mum (bonded), needle size 23 ga, for u 57291-U		63,300	67,300	1315	1315						
89-0810-62	SPME fiber assembly polyacrylate (PA), df 85 mum, for use with autosampler, needle size 23 g 57294-U		66,100	67,300	1315	1315						
89-0810-63	SPME fiber assembly Carboxen/Polydimethylsiloxane (CAR/PDMS), df 85 mum, needle size 23 ga, 57295-U		66,100	67,300	1315	1315						
89-0810-64	SPME fiber assembly Divinylbenzene/Carboxen/Polydimethylsiloxane (DVB/CAR/PDMS), needle size 57298-U		66,100	67,300	1315	1315						
89-0810-81	SPME fiber assembly Polydimethylsiloxane (PDMS), df 100 mum (nonbonded phase), for use with 57341-U		66,100	90,100	1315	1315						
89-0810-85	SPME Fiber Assembly Polydimethylsiloxane/Divinylbenzene (PDMS/DVB), df 65 mum (partially cro 57345-U		66,100	88,400	1315	1315						
89-0810-87	SPME Fiber Holder, for use with CTC CombiPAL, Gerstel MPS2 and Thermo TriPlus Autosamplers 57347-U		64,200	67,800	1315	1315						
89-0810-89	SPME fiber assembly, Carboxen-Polyethylene Glycol (PEG) Coating, needle size 23 ga, df 60 mu 57354-U		72,100	91,000	1315	1315						
89-0843-97	抗生物質-抗真菌剤溶液 (100 x), 細胞生物学用 A5955-100ML A5955		11,800	12,200	1490	1490						
89-0848-02	Accutase (R) solution, sterile-filtered, suitable for cell culture A6964-100ML A6964		9,690	10,100	1490	1490						
89-2368-15	ジフェニルアミン-4-スルホン酸バリウム D4130-25G D4130		50,500	52,300	1500	1500						
89-2368-16	ジフェニルアミン-4-スルホン酸バリウム D4130-50G D4130		98,200	101,000	1500	1500						
89-2368-17	ジフェニルアミン-4-スルホン酸バリウム D4130-5G D4130		16,900	17,500	1500	1500						
89-2372-72	ダルベッコ改変イーグル培地、液体 D6046-500ML D6046		2,780	2,860	1490	1490						
89-2373-46	ダルベッコ改変イーグル培地栄養混合物 F-12/Ham、液体 D6421-500ML D6421		3,520	3,650	1490	1490						
89-2373-48	ダルベッコ改変イーグル培地、液体 D6429-500ML D6429		2,760	2,860	1490	1490						
89-2378-22	ダルベッコリン酸バソファマー、液体 D8537-500ML D8537		2,780	2,860	1490	1490						
89-3089-14	硫酸ゲンタマイシン G3632-250MG G3632		8,580	9,030	1490	1490						
89-3097-37	L-グルタミン、200mM溶液、細胞生物学用 G7513-100ML G7513		5,880	5,690	1490	1490						
89-3099-44	Ala-Gln, 200 mM, solution, sterile-filtered, Biotechnology Performance Certified G8541-100ML G8541		12,100	11,100	1490	1490						
89-3101-68	Geys Balanced Salt Solution, liquid, sterile-filtered, suitable for cell culture G9779-500ML G9779		5,870	6,090	1490	1490						
89-3112-82	HEPES, 1M 溶液、細胞生物学用 H0887-100ML H0887		20,700	21,600	1490	1490						
89-3117-59	HEPES H3375-100G H3375		11,200	11,700	1490	1490						
89-3125-12	ハンクス平衡塩溶液、改変型 H8264-500ML H8264		3,770	3,910	1490	1490						
89-3502-07	イソコフ改変ダルベッコ培地、液体 I3390-500ML I3390		3,230	3,350	1490	1490						
89-3537-91	LB agar, Miller, Microbiologically tested L3027-1KG L3027		53,900	55,800	1500	1500						
89-3537-92	LB agar, Miller, Microbiologically tested L3027-250G L3027		20,200	20,900	1500	1500						
89-3538-84	LBブロス(ミラー)、微生物培養用 L3152-1KG L3152		25,800	26,600	1500	1500</						

品番	商品名	旧価格	新価格	研究用総合 編纂2025 頁	サンクアスト 2025 頁	ナビス2026 頁	サニーフーズ 2024 頁	設備2024 頁	アズビエア 2025 頁	アズツール 2021 頁	サイエンスカ タログ 頁	災害対策カ タログ 頁
89-4809-39	トリトンX-100 X100-100ML X100	11,400	11,900	1490	1490							
89-4809-47	トリトンX-100 X100PC-5ML X100PC	22,200	23,300	1490	1490							
89-4816-43	酵母合成ドロップアウト培地用添加物,酵母培養用 Y0750-20G Y0750	87,900	91,000	1500	1500							
89-4816-63	酵母合成ドロップアウト培地用添加物,酵母培養用 Y1376-20G Y1376	42,000	43,500	1500	1500							
89-4816-67	酵母合成ドロップアウト培地用添加物,酵母培養用 Y1501-20G Y1501	20,000	20,700	1500	1500							
89-4816-76	酵母合成ドロップアウト培地用添加物,酵母培養用 Y1751-20G Y1751	35,700	37,000	1500	1500							
89-4816-78	酵母合成ドロップアウト培地用添加物,酵母培養用 Y1876-20G Y1876	24,700	25,500	1500	1500							
89-4816-80	酵母合成ドロップアウト培地用添加物,酵母培養用 Y2001-20G Y2001	19,400	20,000	1500	1500							
89-6357-38	シリカゲル 13767-500G-R 13767	11,500	11,900	1500	1500							
89-6375-62	クロロホルム 366919-1L 366919	17,100	17,700	1500	1500							
89-6381-27	Water solution, contains 0.1 % (v/v) trifluoroacetic acid, for HPLC 576905-18L 576905	84,100	87,200	1500	1500							
89-6398-90	ウシ胎児血清, 非米国産, 細胞培養用 F7524-500ML	84,300	87,400	1490	1490							
89-6444-50	水酸化カリウム溶液 P4494-50ML P4494	17,300	17,900	1500	1500							
89-6625-19	ウシ胎児血清, オーストラリア産, 細胞培養用 F9423-500ML F9423	226,000	234,000	1490	1490							
89-6630-88	Fetal Bovine Serum, USDA approved, sterile-filtered, suitable for cell culture F0926-500ML F0926	76,800	79,700	1490	1490							
89-6632-02	Fetal Bovine Serum, USA origin, sterile-filtered, suitable for cell culture, suitable for hybr F2442	107,000	111,000	1490	1490							
89-7683-93	Water, for UHPLC, for mass spectrometry 900682-1L 900682	16,800	17,400	1500	1500							
89-7684-42	Supelco 37 Component FAME Mix, certified reference material, TraceCERT(R), in dichloromethane CRM47885	21,800	23,000	1500	1500							
89-7694-22	Water solution, contains 0.1 % (v/v) formic acid, for UHPLC, for mass spectrometry 900687-1L 900687	12,900	13,300	1500	1500							
89-8247-39	Nitinol-core (NIT) SPME Fiber Assembly Polydimethylsiloxane/Divinylbenzene (PDMS/DVB), df 65 57923-U	105,000	111,000	1315	1315							
89-8260-65	Nitinol-core (NIT) SPME Fiber Assembly Carboxen/Polydimethylsiloxane (Carboxen/PDMS), df 75 57907-U	110,000	116,000	1315	1315							
89-8268-54	Nitinol-core (NIT) SPME Fiber Assembly Carboxen/Polydimethylsiloxane (Carboxen/PDMS), for us 57904-U	105,000	111,000	1315	1315							
89-8418-41	30131-U Tenax(R) TA, O.D. x L 1/4 in. x 3 1/2 in., stainless steel TD tube, unconditioned, 35-60 m	136,000	143,000	1315	1315							
9-1191-01	ビタミンシュース(女性用) 22.0cm ホワイトMN410 22W	4,500	5,000			1956						
9-1191-02	ビタミンシュース(女性用) 22.5cm ホワイトMN410 22.5W	4,500	5,000			1956						
9-1191-03	ビタミンシュース(女性用) 23.0cm ホワイトMN410 23W	4,500	5,000			1956						
9-1191-04	ビタミンシュース(女性用) 23.5cm ホワイトMN410 23.5W	4,500	5,000			1956						
9-1191-05	ビタミンシュース(女性用) 24.0cm ホワイトMN410 24W	4,500	5,000			1956						
9-1191-06	ビタミンシュース(女性用) 24.5cm ホワイトMN410 24.5W	4,500	5,000			1956						
9-1191-07	ビタミンシュース(女性用) 25.0cm ホワイトMN410 25W	4,500	5,000			1956						
9-1191-08	ビタミンシュース(女性用) 25.5cm ホワイトMN410 25.5W	4,500	5,000			1956						
9-1192-01	ビタミンシュース(女性用) 22.0cm ピンクMN410 22P	4,500	5,000			1956						
9-1192-02	ビタミンシュース(女性用) 22.5cm ピンクMN410 22.5P	4,500	5,000			1956						
9-1192-03	ビタミンシュース(女性用) 23.0cm ピンクMN410 23P	4,500	5,000			1956						
9-1192-04	ビタミンシュース(女性用) 23.5cm ピンクMN410 23.5P	4,500	5,000			1956						
9-1192-05	ビタミンシュース(女性用) 24.0cm ピンクMN410 24P	4,500	5,000			1956						
9-1192-06	ビタミンシュース(女性用) 24.5cm ピンクMN410 24.5P	4,500	5,000			1956						
9-1192-07	ビタミンシュース(女性用) 25.0cm ピンクMN410 25P	4,500	5,000			1956						
9-1192-08	ビタミンシュース(女性用) 25.5cm ピンクMN410 25.5P	4,500	5,000			1956						
9-1193-01	ビタミンシュース(女性用) 22.0cm パーガンディMN410 22V	4,500	5,000			1956						
9-1193-02	ビタミンシュース(女性用) 22.5cm パーガンディMN410 22.5V	4,500	5,000			1956						
9-1193-03	ビタミンシュース(女性用) 23.0cm パーガンディMN410 23V	4,500	5,000			1956						
9-1193-04	ビタミンシュース(女性用) 23.5cm パーガンディMN410 23.5V	4,500	5,000			1956						
9-1193-05	ビタミンシュース(女性用) 24.0cm パーガンディMN410 24V	4,500	5,000			1956						
9-1193-06	ビタミンシュース(女性用) 24.5cm パーガンディMN410 24.5V	4,500	5,000			1956						
9-1193-07	ビタミンシュース(女性用) 25.0cm パーガンディMN410 25V	4,500	5,000			1956						
9-1193-08	ビタミンシュース(女性用) 25.5cm パーガンディMN410 25.5V	4,500	5,000			1956						
9-1194-01	ビタミンシュース(女性用) 22.0cm ロイヤルブルーMN410 22R	4,500	5,000			1956						
9-1194-02	ビタミンシュース(女性用) 22.5cm ロイヤルブルーMN410 22.5R	4,500	5,000			1956						
9-1194-03	ビタミンシュース(女性用) 23.0cm ロイヤルブルーMN410 23R	4,500	5,000			1956						
9-1194-04	ビタミンシュース(女性用) 23.5cm ロイヤルブルーMN410 23.5R	4,500	5,000			1956						
9-1194-05	ビタミンシュース(女性用) 24.0cm ロイヤルブルーMN410 24R	4,500	5,000			1956						
9-1194-06	ビタミンシュース(女性用) 24.5cm ロイヤルブルーMN410 24.5R	4,500	5,000			1956						
9-1194-07	ビタミンシュース(女性用) 25.0cm ロイヤルブルーMN410 25R	4,500	5,000			1956						
9-1194-08	ビタミンシュース(女性用) 25.5cm ロイヤルブルーMN410 25.5R	4,500	5,000			1956						
7-7356-01	ヘラクレス・ワイヤーカッター 口腔外科、形成外科、手の外科用 OL 601-15-05	123,700	178,500			1310						
7-7356-03	ヘラクレス・ワイヤーカッター 胸膈心臓外科(右手用) OL 603-18	131,500	196,300			1310						
7-7356-04	ヘラクレス・ワイヤーカッター 整形外科用 Φ3.0mm OL 605-23	146,300	218,400			1310						
7-7356-05	ヘラクレス・ワイヤーカッター 整形外科用 Φ3.5mm OL 620-28	209,200	300,300			1310						
7-7356-11	ヘラクレス・ワイヤーカッター 交換用TCブレード(OL 601-15-05用) OL 601-01	42,300	52,700			1310						
7-7356-12	ヘラクレス・ワイヤーカッター 交換用TCブレード(OL 602-18・OL 603-18用) OL 604-01	40,000	54,700			1310						
7-7356-13	ヘラクレス・ワイヤーカッター 交換用TCブレード(OL 605-23用) OL 606-01	40,000	54,700			1310						
7-7356-14	ヘラクレス・ワイヤーカッター 交換用TCブレード(OL 620-28用) OL 607-01	44,800	59,100			1310						
7-4192-11	ネイルマシーン Promed feeling 交換用アタッチメント5本セット 155021	11,300	7,500			646						
7-4192-12	ネイルマシーン Promed feeling 交換用アタッチメント12本セット 198016	14,700	9,800			646						
2-6313-01	プロテティング装置 1704150J1	310,000	320,000	1298	1298							